U.S. DEPARTMENT OF	Office of	
ENERGY	Science	

DOE/SC-ARM-TR-237

Precision Gas System Isotope Analyzer (PGSISO) Instrument Handbook

S Biraud

K Reichl

October 2024



DISCLAIMER

This report was prepared as an account of work sponsored by the U.S. Government. Neither the United States nor any agency thereof, nor any of their employees, makes any warranty, express or implied, or assumes any legal liability or responsibility for the accuracy, completeness, or usefulness of any information, apparatus, product, or process disclosed, or represents that its use would not infringe privately owned rights. Reference herein to any specific commercial product, process, or service by trade name, trademark, manufacturer, or otherwise, does not necessarily constitute or imply its endorsement, recommendation, or favoring by the U.S. Government or any agency thereof. The views and opinions of authors expressed herein do not necessarily state or reflect those of the U.S. Government or any agency thereof.

Precision Gas System Isotope Analyzer (PGSISO) Instrument Handbook

S Biraud K Reichl Both at Lawrence Berkeley National Laboratory

October 2024

How to cite this document:

Biraud, S, and K Reichl. Precision Gas System Isotope Analyzer (PGSISO) Instrument Handbook. 2024. U.S. Department of Energy, Atmospheric Radiation Measurement user facility, Richland, Washington. DOE/SC-ARM-TR-237.

Work supported by the U.S. Department of Energy, Office of Science, Office of Biological and Environmental Research

Acronyms and Abbreviations

Atmospheric Radiation Measurement
American Standard Code for Information Interchange
Bankhead National Forest
Central Calibration Laboratory
U.S. Department of Energy
Data Quality Report
Fourier Transform Infrared
high-resolution transmission molecular absorption line parameters database
input/output
infrared
liters per minute
Multiple Atmospheric Layer Transmission
mercury-cadmium-telluride
National Oceanic and Atmospheric Administration
precise gas system for isotopes
parts per billion
parts per million
quality assurance (data)
quality control/checks (data)
root-mean-square error
Southern Great Plains
Vienna Pee Dee Belemnite
World Meteorological Organization

Contents

Acro	onym	s and Abbreviations	iii
1.0	Instrument Title 1		
2.0	Men	tor Contact Information	. 1
3.0	Insti	ument Design	. 1
	3.1	Instrument/Measurement Theory	. 1
	3.2	Instrument System Functional Diagram	. 3
	3.3	Calibration	. 3
		3.3.1 Corrections	. 4
		3.4.2 Culling Data for Residuals Statistics	. 8
4.0	Defi	nitions and Relevant Information	. 8
	4.1	Primary Variables	.9
	4.2	QC Bit Definitions	.9
	4.3	Calibration Database	.9
5.0	Cita	ble References1	10

Figures

1 Pl	lumbing diagram for PGSISO	;
------	----------------------------	---

Tables

1	Target calibration residual statistics on N Target calibration tank measurements from the SGP	
	deployment AUG 2013 to APR 2023	8
2	WMO CCL scale spans for each species	. 10

1.0 Instrument Title

Precision gas system isotope analyzer (PGSISO)

2.0 Mentor Contact Information

Lead mentor: Sébastien Biraud Lawrence Berkeley National Laboratory Berkeley, California 94720 ph: 510-759-2914 email: <u>SCBiraud@lbl.gov</u>

Associate mentor: Ken Reichl Lawrence Berkeley National Laboratory Berkeley, California 94720 ph: 805-813-1488 email: KReichl@lbl.gov

3.0 Instrument Design

The PGSISO is a combination of several systems in series; (1) gas processing components that pressurize and prevent liquid water entering the sample stream, and (2) the Spectronus FTIR spectrometer, which measures continuous CO_2 , $\delta^{13}C$, $\delta^{18}O$, CH_4 , CO, and N_2O . The Spectronus analyzer has been developed by the Centre for Atmospheric Chemistry at the University of Wollongong, Australia (Griffith et al. 2012), and manufactured by Ecotech Pty Ltd (<u>https://www.ecotech.com/</u>), headquartered in Melbourne, Australia.

Operating parameters internal to the Spectronus are monitored and logged. There are four tower sample inlet channels (for four tower sample heights) and five calibration cylinders that are sequentially selected via Valco multi-port valve. Data are logged by the native Spectronus software. The Spectronus database file is converted to ascii format and further processed to form the pgsiso.b1 datastream.

3.1 Instrument/Measurement Theory

The Spectronus FTIR trace gas analyzer is used to determine simultaneously the concentrations of several trace gases in air or another gas stream by collecting and analyzing the Fourier Transform Infrared (FTIR) spectrum of the gas sample. In the closed-path trace gas analyzer the air is contained in a sample cell that is multi-passed by the infrared (IR) beam. The analyzer consists of separate subsystems controlled by a single computer program:

• The sample handling subsystem, which carries out all sample handling, inlet sample stream selection, cell evacuation, and flushing procedures. This includes a 12-channel digital input/output (I/O) switching capability for switching 12 (or optionally more) solenoid valves, and also an 8-channel analogue-digital converter for logging environmental variables such as temperatures, pressures, and flows.

S Biraud and K Reichl, October 2024, DOE/SC-ARM-TR-237

- The FTIR spectrometer, which produces and records the interferogram and then performs Fourier transformation in order to obtain the infrared spectrum of sampled air. This spectrometer is a Bruker IR cube fitted with a thermoelectrically cooled MCT (mercury-cadmium-telluride) detector (wavenumber range 1500-7800 cm⁻¹). The IR cube is coupled to a 24-meter multi-pass gas cell that contains the air to be analyzed.
- The spectrum analyzer, which provides quantitative analysis of measured spectra immediately after collection, with real-time display.
- Logging of analyzed concentrations, analysis parameters, and analogue input quantities to database file that are parsed to hourly .csv files for ingest by the U.S. Department of Energy (DOE) Atmospheric Radiation Measurement (ARM) user facility.

Optical techniques such as FTIR spectrometry determine the absolute concentrations of individual trace gas species from the quantitative analysis of selected absorption bands or lines in the measured spectrum. The FTIR spectrometer records a broadband mid-IR spectrum that can be analyzed in selected spectral regions for several trace gases. For CH₄, CO, N₂O, and CO₂ combination bands near 3600 cm⁻¹, the low-resolution FTIR spectrometer does not resolve individual isotopic species, and the relevant whole absorption bands are fitted assuming that all isotopologues are present in natural abundance.

For these determinations, natural variations in isotopic abundances do not significantly affect the FTIR measurements. In some cases, such as for ${}^{12}CO_2$ and ${}^{13}CO_2$ analyzed in the v3 bands near 2300 cm⁻¹, individual isotopologue concentrations can be determined independently from the spectrum, allowing the calculation of isotopic ratios.

The quantitative analysis of measured IR absorption spectra uses a computational approach to best-fit regions of the measured spectrum with a calculated spectrum based on knowledge of the sample conditions (pressure, temperature, pathlength), the HITRAN database of spectral line parameters (Rothman et al. 2005, 2009) and a model of the FTIR instrument line shape. Trace gas concentrations and instrument parameters are iteratively adjusted in the calculation to achieve best fit. This approach uses the program MALT (Multiple Atmospheric Layer Transmission; Griffith 1996) to calculate spectra and is described in more detail in Griffith et al. 2012 and the MALT user's manual.

From each measured spectrum, MALT retrieves the concentration of each target gas and converts to a mole-fraction scale using the measured total pressure and temperature of the sample. The raw FTIR scale depends on the accuracy of HITRAN parameters, the MALT spectrum model, the optical pathlength, temperature, pressure, and other parameters. It is usually within a few percent of established standard reference scales (Smith 2011).

3.2 Instrument System Functional Diagram



Figure 1. Plumbing diagram for PGSISO.

3.3 Calibration

Measurements of species are calibrated by measuring four span cylinders about every seven days, with measurements of one Target cylinder at the same pace, and another once per day. Linearly interpolated gain and offset coefficients for each species are applied to the sample data that undergo several stages of corrections prior to calibration. Notable regimes for reference calibration are:

- August 2013 deployment to September 2015 3-point calibration, linear interpolation.
- September 2015 to April 2023 4-point calibration, linear interpolation.
- January 2024 to present 5-point calibration, quadratic interpolation.

Below are conversions, corrections, and calibration steps included in the a1-level to b1-level ingest processing.

3.3.1 Corrections

H₂O Correction

Water vapor displacement correction must occur for all species measured. The H₂O correction goes as:

$$\chi' = \frac{\chi}{\left(1 - \frac{q}{10^6}\right)}$$

where χ is the concentration of current sample measurement of species ¹²CO₂, ¹³CO₂, ¹⁸O¹²C¹⁶O, CH₄, CO, N₂O. and q is the concentration of H₂O at the coincident to measurement of χ .

Residual Correction

How and whether the residual correction is applied are determined based on measurement modes operated throughout the key time periods pertaining to the ARM Southern Great Plains (SGP) observatory deployment:

- August 2013 deployment to September 2015
 - Tower and Target tank (inlet 1-4, and 10, respectively)
 - $\circ \quad P_{resid} = P_{evac2} + 83.2 + 8.3 + 4.0$
 - where P_{evac2} is the evacuated pressure of previous sample measurement
 - Cal tanks (inlet 5-9)
 - \circ P_{resid} = 4.0
- September 2015 to January 2017
 - Tower (inlet 1-4)
 - $\circ \quad P_{\text{resid}} = P_{\text{evac2}} + 8.3 + 4.0$
 - Target and Cal tanks (inlet 5-9, 10, 16)

$$\circ$$
 P_{resid} = 4.0

• January 2017 to present \rightarrow No residual correction.

Data collected for the upcoming Bankhead National Forest (BNF) deployment in Alabama also do not require a residual correction.

Correction for the presence of residual air sample from previous measurement is required, depending on whether the measurement cell is flushed or not flushed with current sample air prior to measurement of the current sample.

$$\chi' = \frac{\left(\chi P_{cell} - \chi_{prev} P_{resid}\right)}{\left(P_{cell} - P_{resid}\right)}$$

where χ is the H₂O corrected concentration of sample measurement of species ¹²CO₂, ¹³CO₂, ¹⁸O¹²C¹⁶O, CH₄, CO, N₂O; χ_{prev} is the H₂O corrected concentration from previous adjacent sample measurement; P_{resid} is the residual pressure of sample from previous measurement; P_{cell} is the cell pressure of current sample measurement.

Cross-Sensitivity Correction

This stage is for converting data to standard conditions. To first order, water vapor interference and sample pressure and temperature effects on the spectrum are included in the MALT calculation and spectrum fit. However, due to small inaccuracies and uncertainties in the measured spectrum, HITRAN line parameters, the spectral model, and the measured temperature and pressure, small dependences of apparent (raw) measured mole fractions on these parameters may remain.

$$\chi' = \chi - \frac{d\chi}{dQ}(Q - Q_0) - \frac{d\chi}{dP}(P - P_0) - \frac{d\chi}{dT}(T - T_0) - \frac{d\chi}{dF}(F - F_0) - \frac{d\chi}{dC}(C - C_0)$$

where χ is the H₂O corrected and if applicable, residual corrected concentration of sample measurement of species ¹²CO₂, ¹³CO₂, ¹⁸O¹²C¹⁶O, CH₄, CO, N₂O; d χ d* is the cross-sensitivity coefficient for the species; Q is the concentration of H₂O and Q₀ = 0 ppm; P is the cell pressure and P₀ = 1100 mbar; T is the cell temperature and T₀ = 35°C; F is the cavity sample flow and F₀ = 1 LPM; C is the concentration of CO₂ and C₀ = 400 ppm.

Rescale CO₂ Isotope Measurements from HITRAN to VPDB scale

To allow for calculation of δ^{13} C on the Vienna Pee Dee Belemnite (VPDB) scale:

$$\label{eq:constraint} \begin{split} ^{12} \chi_{vpdb} &= ^{12} \chi' \ x \ 1.000150 \\ ^{13} \chi_{vpdb} &= ^{13} \chi' \ x \ 1.005280 \\ ^{18} \chi_{vpdb} &= ^{18} \chi' \ x \ 0.96319 \end{split}$$

where ${}^{12}\chi'$, ${}^{13}\chi'$ and ${}^{18}\chi'$ are the H₂O-corrected, residual-corrected, cross-sensitivity-corrected concentration of scaled mole fractions for ${}^{12}CO_2$, ${}^{13}CO_2$, and ${}^{18}O^{12}C^{16}O$ calculated from the HITRAN scale, respectively.

3.3.2 Reference Calibration

Calibration coefficients are determined by linear or quadratic regression by measuring known values of species for a set of calibration tanks. For linear calibration, gain and offset coefficients are calculated for each species. For quadratic calibration, curve, gain, and offset coefficients are calculated for each species.

The linear calibration for data from Aug 2013 to Apr 2023 goes as:

$$\chi_{\rm corr} = (\chi' - b) / a$$

where χ' is the measured concentration of each species ${}^{12}CO_2$, ${}^{13}CO_2$, ${}^{18}O^{12}C^{16}O$, CH₄, CO, N₂O with relevant corrections from Section 3.3.1 applied, *b* is the offset coefficient, *a* is the gain coefficient, and χ_{corr} is the calibrated concentration of each species.

The quadratic calibration for data from Jan 2024 to present goes as:

$$\chi_{\rm corr} = c \, \chi'^2 + a \, \chi' + b$$

where χ' is the measured concentration of each species ${}^{12}CO_2$, ${}^{13}CO_2$, ${}^{18}O^{12}C^{16}O$, CH₄, CO, N₂O with relevant corrections from Section 3.3.1 applied, *c* is the curve coefficient, *b* is the offset coefficient, *a* is the gain coefficient, and χ_{corr} is the calibrated concentration of each species.

CO2 Isotope Calibration

Known values of calibration cylinders are reported for δ^{13} C, δ^{18} O, and CO₂, while scaled mole-fraction values are reported by the analyzer as χ'_{626} , χ'_{636} , χ'_{628} , for the ¹⁶O¹²C¹⁶O, ¹⁶O¹³C¹⁶O, and ¹⁸O¹²C¹⁶O isotopologues respectively, and must first be converted. The steps to calculate calibrated CO₂, δ^{13} C, and δ^{18} O are:

1. Convert the reported known values of CO₂, δ^{13} C, and δ^{18} O into the scaled mole-fraction values for χ_{626} , χ_{636} and χ_{628} as:

$$\chi_{626} = \frac{y_{c02}}{R_{sum} x_{626}}$$
$$\chi_{636} = \frac{(1 + \delta^{13}) y_{c02}}{R_{sum} x_{626}}$$
$$\chi_{628} = \frac{(1 + \delta^{18}) y_{c02}}{R_{sum} x_{626}}$$

where y_{co2} is the known value of the mole fraction of CO₂, δ^{13} is shorthand for the known value of the isotope ratio δ^{13} C/1000, similarly for δ^{18} as the known value of the isotope ratio δ^{18} O/1000, R_{sum} is the approximated sum of the dominant CO₂ isotopologues as:

$$R_{sum} = (1 + {}^{13}r)(1 + {}^{17}r + {}^{18}r)^2$$

where,

$${}^{13}r = {}^{13}r_{ref}(1 + \delta^{13})$$

$${}^{18}r = {}^{18}r_{ref}(1 + \delta^{18})$$

$${}^{17}r = {}^{17}r_{ref}(1 + (0.528 \ \delta^{18}))$$

and,

 ${}^{13}r_{ref} = 0.111802$, the abundance ratio of ${}^{16}O{}^{13}C{}^{16}O$ to ${}^{16}O{}^{12}C{}^{16}O$, ${}^{18}r_{ref} = 0.00208835$, the abundance ratio of ${}^{18}O{}^{12}C{}^{16}O$ to ${}^{16}O{}^{12}C{}^{16}O$, ${}^{17}r_{ref} = 0.0003931$, the abundance ratio of ${}^{17}O{}^{12}C{}^{16}O$ to ${}^{16}O{}^{12}C{}^{16}O$,

and,

 $x_{626} = 0.984054$, the abundance of ¹⁶O¹²C¹⁶O in reference to VPDB

- 2. Calculate the calibration coefficients independently for χ_{626} , χ_{636} and χ_{628} in the same manner as the other species.
- 3. Calculate calibrated CO₂, δ^{13} C and δ^{18} O using χ_{corr} , the calibrated scale mole-fraction values from each isotopologue:

$$\delta^{13}C = \left(\frac{\chi'_{636,corr}}{\chi'_{626,corr}} - 1\right) \times 1000\%$$

$$\delta^{18}O = \left(\frac{\chi'_{628,corr}}{\chi'_{626,corr}} - 1\right) \times 1000\%$$

$$y_{co2} = \chi'_{626,corr}R_{sum}x_{626}$$

where, $\chi'_{626,corr}$, $\chi'_{636,corr}$ and $\chi'_{628,corr}$ are the calibrated scaled mole fractions of the isotopologues ${}^{16}O^{12}C^{16}O$ and ${}^{16}O^{13}C^{16}O$ and ${}^{18}O^{12}C^{16}O$ respectively, with relevant corrections from Section 3.3.1 applied, y_{co2} is the calibrated total isotopic mole fraction of CO₂, and

$$R_{sum} = (1+{}^{13}r)(1+{}^{17}r+{}^{18}r)^2$$

where,

$${}^{13}r = {}^{13}r_{ref} (1 + \delta^{13}C/1000)$$

$${}^{18}r = {}^{18}r_{ref} (1 + \delta^{18}O/1000)$$

$${}^{17}r = {}^{17}r_{ref} (1 + (0.528 \ \delta^{18}/1000))$$

where $\delta^{13}C$ and $\delta^{18}O$ are the calibrated isotope ratios

3.4 Uncertainty

Overall instrument uncertainty is characterized by the variability in the calibration residual from measurements of Target cylinders compared to the known secondary values of each species, across the record from October 2013 to present (as of the date of this handbook).

This allows for the uncertainty to encompass components of uncertainty that are caused by systematic and/or operational errors that may or may not be characterized and that come from the actual operation of the system.

3.4.1 Uncertainty from Calibration Residuals

The Target calibration residual is defined as: Residual = (Calibrated, Measured) – (Known Value).

The table below represents uncertainty metrics used by mentors to assess instrument performance, where columns are defined as:

• MEAN = mean of difference measured and known quantities as $\chi_{resid} = \chi_{meas} - \chi_{known}$

where χ are subscripted as measured calibrated or known values of species CO₂, δ^{13} C, δ^{18} O, CH₄, CO, and N₂O

- STDERR = the standard error of the mean of residuals
- RMSE = the root-mean-square error calculated as,

 $\sqrt{(\sum \chi_{\text{resid}}^2 / N)}$

where the summation occurs over the square of residual from of N Target tank measurements used in the assessment.

S Biraud and K Reichl, October 2024, DOE/SC-ARM-TR-237

	MEAN	STDERR	RMSE	Ν
CO ₂ (ppm)	-0.1	0.0	0.2	1313
δ ¹³ C (‰)	0.1	0.0	0.2	1350
δ ¹⁸ O (‰)	-0.1	0.0	0.5	1404
CH ₄ (ppb)	-0.1	0.0	0.3	1329
CO (ppb)	1.3	0.0	1.6	1325
N ₂ O (ppb)	0.1	0.0	0.2	712

Table 1.Target calibration residual statistics on N Target calibration tank measurements from the SGP
deployment AUG 2013 to APR 2023.

3.4.2 Culling Data for Residuals Statistics

Not all Target cylinder measurements are used in the assessment of instrument uncertainty. In addition to the diagnostic QC tests the measurements undergo (Section 4.2), specified periods of measurements of the Target tank are excluded from statistics using machine-readable documentation describing possible or confirmed periods of instrument instability coincident with those measurements. This machine-readable documentation is a mentor-compiled .csv file of all of the Data Quality Reports (DQRs) associated with the datastream and is used as a configuration file in the data ingest process. Examples of criteria for manually identified periods of instrument instability include:

- If H₂O measurements are abnormally high or variable
- Calibration coefficient outliers or periods of systematic deviation from expected values; significant outliers suggest the instrument was not in optimal condition and may be coincident with the previous criteria.

For the case of coefficients, when data is being reprocessed and coefficients are recalculated, the QC state of coefficients corresponding to the DQR entries are overwritten as "missing value" to then be excluded from linear interpolation of coefficients.

4.0 Definitions and Relevant Information

The b1-level datastreams are: sgppgsisoC1.b1 bnfpgsisoS10.b1

Auxiliary and input datastreams associated with b1-level datastreams;

- for SGP:
 - sgppgsisoC1.a1
 - sgppgsisoC1.b0
 - sgppgsisocoeffC1.b1
- for BNF:
 - bnfpgsisoS10.a1
 - bnfpgsisoS10.b0
 - bnfpgsisocoeffS10.b1

4.1 Primary Variables

At b1-level, primary variable names are:

- CO2_AVG_CORR
- CH4_AVG_CORR
- CO_AVG_CORR
- N2O_AVG_CORR
- DELTA13C
- DELTA180

Associated QC variables:

- CO2_AVG_CORR_QC
- CH4_AVG_CORR_QC
- CO_AVG_CORR_QC
- N2O_AVG_CORR_QC
- DELTA13C_QC
- DELTA180_QC

Note:

- The "best data" are those with *_QC values of 0 for tower sample or Target measurements
- Tower sample data correspond to VALCOPOSITION variable equal to integers 1, 2, 3, or 4, for the four tower sample heights respectively.

4.2 QC Bit Definitions

```
flag_method = "bit";
bit_1_description = "Value is equal to missing_value -9999."; 1
bit_2_description = "Value is less than the valid_min."; 2
bit_3_description = "Value is greater than the valid_max."; 4
bit_4_description = "Value for *_slope is missing value"; 8
bit_5_description = "Value for *_slope_err is missing value"; 16
bit_6_description = "Value for *_slope is below valid_min."; 32
bit_7_description = "Value for *_slope is above valid_max."; 64
bit_8_description = "Value for *_slope is significant with respect to a slope of zero, where
(abs(co2_slope-co2_slope_err) > 0)"; 128
bit_9_description = "QC state for at least one of the mentor-defined auxiliary variables is not
equal to 0"; 256
bit_10_description = "Measurement follows a tank for which residual correction is not a
sufficient correction"; 512
bit_11_description = "Measurement of N2 Purge Tank"; 1024
```

4.3 Calibration Database

Calibration tank values are transcribed to a machine-readable configuration file used in ARM ingest processing. The file is updated by the mentor each time the calibration tanks are changed out – about once

S Biraud and K Reichl, October 2024, DOE/SC-ARM-TR-237

every five years. The values are provided by the World Meteorological Organization (WMO)'s Central Calibration Laboratory (CCL) at the U.S. National Oceanic and Atmospheric Administration (NOAA)'s Earth System Research Laboratory. The calibration scales used for each of the species for current installed calibration cylinders (July 2018) are given below in Table 2.

	Calibration Scale
CO ₂	WMO-CO ₂ -X2019
CH ₄	WMO-CH ₄ -X2004A
СО	WMO-CO-X2014A
N ₂ O	WMO-N ₂ O-X2006A

Table 2.WMO CCL scale spans for each species.

Contact the mentor for access to calibration certificates of each species.

5.0 Citable References

Griffith, DWT. 1996. "Synthetic calibration and quantitative analysis of gas-phase FT-IR spectra." *Applied Spectroscopy* 50(1): 59–70, <u>https://doi.org/10.1366/0003702963906627</u>

Griffith, DWT, NM Deutscher, G Kettlewell, M Riggenbach, C Caldow, and S Hammer. 2012. "A Fourier transform Infrared trace gas analyser for atmospheric applications." *Atmospheric Measurement Techniques* 5(10): 2481–2498, https://doi.org/10.5194/amt-5-2481-2012

Griffith, DWT. 2018. "Calibration of isotopologue-specific optical trace gas analyzers: a practical guide." *Atmospheric Measurement Techniques* 11(11): 6189–6201, <u>https://doi.org/10.5194/amt-11-6189-2018</u>

Rothman, LS, D Jacquemart, A Barbe, DC Benner, M Birk, LR Brown, MR Carleer, C Chackerian, K Chance, V Dana, VM Devi, J-M Flaud, RR Gamache, A Goldman, J-M Hartmann, KW Jucks, AG Maki, J-Y Mandin, ST Massie, J Orphali, A Perrin, CP Rinsland, MAH Smith, J Tennyson, RN Tolchenov, RA Toth, JV Auwera, P Varanasi, and G Wagner. 2005. "The HITRAN 2004 molecular spectroscopic database." *Journal of Quantitative Spectroscopy and Radiative Transfer* 96(2): 139–204, https://doi.org/10.1016/j.jqsrt.2004.10.008

Rothman, LS, IE Gordon, A Barbe, DC Benner, PE Bernath, M Birk, V Boudon, LR Brown, A Campargue, JP Champion, K Chance, LH Coudert, V Dana, VM Devi, S Fally, JM Flaud, RR Gamache, A Goldman, D Jacquemart, I Kleiner, N Lacome, WJ Lafferty, JY Mandin, ST Massie, SN Mikhailenko, CE Miller, N Moazzen-Ahmadi, OV Naumenko, AV Nikitin, J Orphal, VI Perevalov, A Perrin, A Predoi-Cross, CP Rinsland, M Rotger, M Simeckova, MAH Smith, K Sung, SA Tashkun, J Tennyson, RA Toth, AC Vandaele, and J Vander Auwera. 2009. "The HITRAN 2008 molecular spectroscopic database." *Journal of Quantitative Spectroscopy and Radiative Transfer* 110(9-10): 533–572, https://doi.org/10.1016/j.jqsrt.2009.02.013 Smith, TEL, MJ Wooster, M Tattaris, and DWT Griffith. 2011. "Absolute accuracy and sensitivity analysis of OP-FTIR retrievals of CO₂, CH₄ and CO over concentrations representative of 'clean air' and 'polluted plumes'." *Atmospheric Measurement Techniques* 4(1): 97–116, <u>https://doi.org/10.5194/amt-4-97-2011</u>

World Meteorological Organization (WMO). 2014. "GAW Report No. 213." From 17th WMO/IAEA Meeting of Experts on Carbon Dioxide, other Greenhouse Gases and Related Tracers Measurements Techniques. Beijing, China.

http://www.wmo.int/pages/prog/arep/gaw/documents/Final_GAW_213_web.pdf



www.arm.gov



Office of Science