AN AIRBORNE DIGITAL HOLOGRAPHIC INSTRUMENT FOR MEASURING THE SPATIAL DISTRIBUTION AND LOCAL SIZE DISTRIBUTIONS OF CLOUD PARTICLES: HOLOGRAPHIC DETECTOR FOR CLOUDS 2 (HOLODEC 2)

> JACOB P. FUGAL, SCOTT SPULER EARTH OBSERVING LABORATORY NCAR, BOULDER, CO USA & RAYMOND A. SHAW PHYSICS DEPARTMENT, MICHIGAN TECH HOUGHTON, MI USA



HOLODEC (Holographic Detector for Clouds) is an airborne instrument that measures the size, shape, and relative 3D position of cloud particles using digital inline holography.

SCIENCE QUESTIONS FOR HOLODEC 2

- How do local cloud particle size distributions vary inside cloud regions (edge, top, core, base), by cloud age, cloud type?
- How are cloud particles spatially distributed on sub-cm scales due to mixing, entrainment, and turbulent processes?
- How are ice and liquid water particles spatially distributed or partitioned in mixed-phase regions of cloud?

HOLOGRAPHIC MEASUREMENT UNIQUENESS AND ADVANTAGES

- Sample volume is the shape and size of a finger as opposed to a order 100 m long ribbon or pencil shape. Allows localized measurements.
- Sample volume size is insensitive to a particle size dependent depth of focus or air speed. More accurate particle concentration estimates.
- Measures three-dimensional position of particles. Allows view of cloud texture and structure and more accurate rejection of shattered particles.

MEASUREMENT PRINCIPLE: DIGITAL IN-LINE HOLOGRAPHY



Hologram (simulated) 1/20th actual size

HOLOGRAM RECONSTRUCTED

HOLODEC 1: PARAMETERS

- Minimum particle size: ~22 μm
- Sample volume size: 4.8 x 3.6 x 30 mm or ~ 0.5 cm³
- Camera frame rate: ~15 Hz
- Volume sample rate:
 ~7 cm³ s⁻¹
- Hologram processing time: ~22 min

HOLODEC 1

RECONSTRUCTED ICE CRYSTALS

DETECTING ICE PARTICLE SHATTERING



DETECTING ICE PARTICLE SHATTERING



NUMBER DENSITIES





SIZE DISTRIBUTIONS

CURRENT HOLODEC 2 PARAMETERS



- Minimum particle size: ~5 μm diameter
- Sample volume: ~ 10 cm³
- Sample Volume Spacing @ 200 m/s: 66 m
- Volume sample rate: ~30 cm³ s⁻¹
- Hologram processing time: ~ 2 hours

HOLODEC 2 PARAMETERS WITH A FAST CAMERA



- Minimum particle size: ~5 μm diameter
- Sample volume: ~ 1 cm³
- Sample Volume Spacing @ 200 m/s: 40 cm
- Volume sample rate: ~ 400 cm³ s⁻¹

HOLODEC 2 PARAMETERS IN GENERAL

- Size uncertainty is order 1 pixel width or 3 µm. Can be improved with better hologram processing algorithms.
- Particle phase determination is by shape so works for particles order 10 pixel widths or 30 µm diameter.
- Low power, lightweight, standard PMS canister form factor except for Gig-E connection.
- Autonomous operation possible. Instrument is essentially a airborne digital camcorder.

CURRENT PLANS

- Planned test flights in Summer 2009
- Intended hologram processing computer: US NSF sponsored Teragrid supercomputer Lincoln, heterogenous (CPU & GPU) compute cluster.

SUMMARY

- Digital in-line holography allows in-situ measurements of three-dimensional position, and two-dimensional profile and size of each cloud particle in a prism-shaped well-defined sample volume.
- Using these measurements HOLODEC 2 will be able to measure local size distributions of cloud particles without shattering artifacts, or sample volume size and air speed uncertainties.
- We will also be able to examine the fine scale structure of clouds to learn more about mixing, entrainment, and turbulence processes (particle clustering), and the fine-structure of mixed-phase clouds.